### NCH - PEDMATCH

Integration Overview, v1.0-draft

Last revised, 5/9/2016, ESK


## Terminology

- Collection Grouping ID ("CGID")
  - Shared by specimens of a given type taken from a patient in one encounter/surgery/draw/etc.
  - Same CGID means specimens are "equivalent" from a PED-MATCH standpoint
  - Specimen messages sent to MATCHBOX will include the CGID
  - CGIDs will be "inherited" by the various aliquots/derivatives from the original specimens
  - Specimens obtained through follow-up requests by BPC (e.g., inadequate/insufficient received) will keep the original CGID if specimens are known to come from the same original encounter/surgery/block/etc.

## Collection Grouping Examples

### Scenario 1

- First two specimens grouped together because attributes match ("equivalent" from PED-MATCH standpoint); New (shared) CGID
- Third specimen was collected separately from first two; New CGID
- Fourth specimen differs by type; New CGID



time

# Collection Grouping Examples

#### Scenario 2

- First two specimens differ in collection date; Each receives a unique CGID
- Third specimen—though arriving at the BPC days later—reuses the first specimen's CGID since they are "equivalent" from a PED-MATCH standpoint
  - Late-arriving specimens (e.g., delays in shipping/accessioning)
  - Cases where BPC requests additional portion of existing specimen (e.g., inadequacy) and site is able to indeed provide more from that original specimen



# Collection Grouping Examples

#### Scenario 3

- First two specimens differ in collection date; Each receives a unique CGID
- Third specimen arrives at the BPC days later but does NOT reuse the first specimen's
   CGID since they are not from the same collection and thus not "equivalent" from a PED-MATCH standpoint
  - Cases where BPC requests additional portion of existing specimen (e.g., inadequacy) but site is unable to provide more of the original specimen, instead substituting from another specimen/draw/etc.




### "Specimen Received" Integration Message

- DIRECTION: Outbound from BPC to MATCHBOX
- <u>TIMING</u>: Sent when new specimens of interest are detected in STARS
  - Specimen is accessioned under APEC1621
  - Patient is enrolled in APEC1621
  - Specimen type of Blood or Tumor Tissue
  - Accessioning paperwork QC is successful

### CONTENT:

- Technical Message Headers (msg\_guid, msg\_dttm)
- List of specimens
  - PED-MATCH attributes (patient\_id, collection\_grouping\_id, type, disease\_status, study\_id, received\_ts, collection\_ts)
  - Additional (troubleshooting) attributes (stars\_patient\_id, stars\_specimen\_id, stars\_specimen\_type, qc\_ts)
- Timestamps conveyed in GMT/ISO-8601

#### BPC / PED-MATCH - Process and Integration Overview - Specimens Received (Success)



"Specimen Received" Integration Message

```
1 - {
2 +
        "msg guid": "5c64192f-8a25-4874-9db6-fd55c398822d",
4
        "msg dttm": "2016-04-25T18:42:13+00:00"
5
                                                                                                             4
 6 +
      "specimens_received": [
7 +
8
           "patient id": "17364",
9
           "collection_grouping_id": "CGID-123456",
10
           "type": "Blood",
11
          "disease_status": "Normal",
12
          "collection_ts": "2016-04-25T14:17:11+00:00",
13
          "received ts": "2016-04-25T15:17:11+00:00",
14
          "study_id": "ABC1234",
15 -
          "internal_use_only": {
16
            "stars patient id": "ABCXYZ",
17
             "stars specimen id": "ABCXYZ-ABC123",
            "stars specimen type": "Fresh Blood",
18
19
             "ac ts": "2016-04-25T16:21:34+00:00"
20
21
22 +
23
           "patient id": "23123",
24
          "collection grouping id": "CGID-77864",
           "type": "Tissue".
25
           "disease status": "Diseased",
26
           "collection ts": "2016-04-25T14:17:11+00:00",
27
           "received ts": "2016-04-25T15:17:11+00:00",
28
           "study id": "ABC1234".
29
30 +
           "internal use only": {
            "stars patient id": "QWERTY",
31
32
             "stars specimen id": "QWERTY-XYZ222",
33
             "stars specimen type": "Fresh Tissue Primary",
34
             "qc ts": "2016-04-25T16:21:34+00:00"
35
36
37 +
38
           "patient_id": "23123",
           "collection_grouping_id": "CGID-77864",
39
40
           "type": "Tissue",
41
           "disease status": "Diseased",
42
           "collection ts": "2016-04-25T15:17:11+00:00",
43
           "received_ts": "2016-04-25T16:17:11+00:00",
44
           "study_id": "ABC1234";
45 +
           "internal_use_only": {
46
            "stars_patient_id": "QWERTY",
47
            "stars_specimen_id": "QWERTY-XYZ223",
48
            "stars_specimen_type": "Fresh Tissue Primary",
            "qc_ts": "2016-04-25T17:21:35+00:00"
49
50
51
52 +
           "patient_id": "23123",
53
54
           "collection grouping id": "CGID-77865",
55
           "type": "Blood",
56
           "disease status": "Normal",
57
           "collection ts": "2016-04-25T14:17:11+00:00",
58
          "received ts": "2016-04-25T15:17:11+00:00",
          "study id": "ABC1234",
59
          "internal_use_only":
60 +
```

```
    □ ▼ object {2}

      ▼ header {2}
             msg guid: 5c64192f-8a25-4874-9db6-fd55c398822d
             msg dttm: 2016-04-25T18:42:13+00:00
      ▼ specimens received [4]
                patient id: 17364
                collection grouping id: CGID-123456
                type : Blood
                disease status : Normal
                collection_ts: 2016-04-25T14:17:11+00:00
                received ts: 2016-04-25T15:17:11+00:00
                study id: ABC1234
             ▼ internal use only {4}
                   stars patient id : ABCXYZ
                   stars_specimen_id : ABCXYZ-ABC123
                   stars_specimen_type : Fresh Blood
                   qc ts: 2016-04-25T16:21:34+00:00
          ▶ 1 {8}
          ▼ 3 {8}
                patient_id : 23123
                collection_grouping_id : CGID-77865
                type : Blood
                disease status : Normal
                collection ts: 2016-04-25T14:17:11+00:00
                received ts: 2016-04-25T15:17:11+00:00
                study id: ABC1234
             ▼ internal use only {4}
                   stars patient id : QWERTY
                   stars specimen id: QWERTY-XYZ333
                   stars_specimen_type : Fresh Blood
                   gc ts: 2016-04-25T16:21:34+00:00
```

### To be continued...

- Definition of specific attributes for grouping specimens for CGID derivation
- Semantics and messaging scenarios for "specimen failed" messaging

## GitHub Repo

- Public repository
- Read-only access to NCH design documents, sample messages, etc.
- Supports versions and branching
- https://github.com/NCH-BPC-Informatics/pedmatcharch-docs

